

When the cell supernatant of ATL-72b(-) is added to Vero cells which are indicator cells for mycoplasma tests, <u>unidentified autonomous proliferative particles are observed after several days.</u> (*The particles are not observed in the culture of ATL-72b(-) cells alone.) With our mycoplasma tests (DNA staining test and PCR test) and Bacteria/Fungi PCR test, they all indicated negative results.





*Similar phenomena have been reported with other cell lines of other cell banks.

Harasawa, R. et al., In vitro Cell. Dev. Biol. Anim. 42: 13-15, 2006.

At present, it is unknown why these particles increase, and there is no evidence that these are live particles.