

## **Media preparation for BHK/T7-9 cells (RCB4942)**

MEM + 5% FBS + 10% Tryptose Phosphate Broth

Ex) Preparation of culture media : 100mL

85mL MEM

5mL FBS

10mL Tryptose Phosphate Broth\*

\*Preparation of Tryptose Phosphate Broth (BD, Difco 260300) : 1L

1. Weigh 29.5g of Tryptose Phosphate Broth powder and measure 1L of purified water.
2. Dissolve the powder in the water by swirling.
3. Autoclave at 121°C for 15 minutes.