Application consideration

RCB5005: FKH-1

Please keep in mind the following points:

- · The cells will grow slowly.
- Don't change the medium on the next day of thawing.
 Keep observing cells for 1 week without medium change.
- Since the activity of the cytokine to be added varies by the manufacturer, please be aware of the additive concentration.
 If the cells are not growing well, please prepare and observe several dishes with different concentrations of cytokine.
- *Cytokines are consumed over time.

When adding cytokines to the medium, please calculate the required amount based on the total volume of the medium in the culture vessel.

【Technical support】 e-mail : cellqa@brc.riken.jp