

Month/Day/Year

To whom it may concern,

The following cell line was provided to the user Dr.?????, on Month/Day/Year:
Cell Name: ????? (RCB ?????)

Prior to the provision, the following validation of the above cell line was performed:

DNA barcoding analysis:

DNA barcoding analysis demonstrated that this cell line was derived from ????.

Mycoplasma infection:

Vero cells were cultured for 5 days in the presence of the culture supernatant of the cell line. The cultured Vero cells were stained with Hoechst 33258 for DNA staining, and no extra-nuclear DNA staining was detected indicating that the presence of mycoplasma in the culture supernatant of the cell line was below the detection level.

(sign)

Kazuaki NAKAMURA, Ph.D.

Director of the Cell Engineering Division RIKEN BioResource Research Center
Koyadai 3-1-1, Tsukuba, Ibaraki, 305-0074, Japan

URL: <https://cell.brc.riken.jp/en/>

E-mail: cellbank@brc.riken.jp

FAX: +81-29-836-9130