

## Application consideration

### RCB3010 : RPMI8226

Please keep in mind the following points :

- This cell line grows as a combination of adherent and floating cells.
- The cells typically require approximately one week to proliferate after seeding.  
(Please note that the proliferation rate may vary depending on the type of serum used.)
- Basically, the cells should be cultured while observing their condition.
- If the culture medium shows signs of discoloration (e.g., orange or yellow) three to four days after the frozen stock is seeded, please add fresh medium equal to about half of the total volume.
- During passaging, please collect both suspension and adherent cells.  
(Make sure to confirm detachment under a microscope.)

\*For details on how to culture the cells after delivery, please refer to our website below.  
[https://cell.brc.riken.jp/en/wp-content/uploads/pdf/RIKEN BRC CELL BANK Information Sheet.pdf](https://cell.brc.riken.jp/en/wp-content/uploads/pdf/RIKEN_BRC_CELL_BANK_Information_Sheet.pdf)

